



SPECIMEN TYPE	WOLFE Test	TESTING FOR	WINDOWS OF DETECTION	COLLECTION CONSIDERATIONS	TIME REQUIRED FOR COLLECTIONS	RESULTS AVAILABLE	OTHER IMPORTANT CONSIDERATIONS
INSTANT URINE TEST	<b>RealityCHECK Premium Urine Test Cup</b>	3-12 drugs of abuse	<ol style="list-style-type: none"> <li>As early as two hours after use.</li> <li>On average up to 72 hours after use.</li> <li>Depends upon individual metabolism and drug.</li> </ol>	Collection site or on-site bathroom.	<p>If on-site, from 15 minutes up to three hours.</p> <p>Off-site, depends upon distance traveled.</p>	Five minutes	<ol style="list-style-type: none"> <li>Two-step process activated by collector to ensure tamper-proof collection.</li> <li>Zero exposure to the specimen.</li> <li>Has the most features.</li> </ol>
INSTANT URINE TEST	<b>RealityCHECK One-Step Drug Screen Card</b>	1-12 drugs of abuse	<ol style="list-style-type: none"> <li>As early as 2 hours after use.</li> <li>On average up to 72 hours after use.</li> <li>Depends on individual metabolism and drug.</li> </ol>	Collection site or on-site bathroom.	<p>If on-site, from 15 minutes up to three hours.</p> <p>Off-site, depends upon distance traveled.</p>	Five minutes	<ol style="list-style-type: none"> <li>Least expensive alternative.</li> </ol>
LAB-BASED ORAL FLUID TEST	<b>RealityCHECK Premium Oral Fluid Collector</b>	Seven drugs of abuse including THC parent compound & Ecstasy	<ol style="list-style-type: none"> <li>As early as right away.</li> <li>On average up to 36 hours after use.</li> <li>Depends upon drug.</li> </ol>	<p>On-site, no restrictions.</p> <p>Avoids gender-specific collection issues.</p>	If on-site, from 5 to 15 minutes.	As early as same day received by lab.	<ol style="list-style-type: none"> <li>Neat specimen improves accuracy of results.</li> <li>Helps eliminate adulteration concerns.</li> <li>Virtually eliminates quantity not sufficient (QNS) rejections.</li> <li>Eliminates personnel issues associated with HR or supervisors testing employees.</li> </ol>
LAB-BASED URINE TEST	<b>Urine Lab Testing</b>	1-12 drugs of abuse	<ol style="list-style-type: none"> <li>As early as 2 hours after use.</li> <li>Up to 72 hours after use.</li> <li>Depends on individual metabolism and drug.</li> </ol>	Collection site or on-site bathroom.	<p>If on-site, from 15 minutes up to three hours.</p> <p>Off-site, depends on distance traveled.</p>	As early as same day received by lab.	<ol style="list-style-type: none"> <li>Eliminates personnel issues associated with HR or supervisors testing employees.</li> </ol>
LAB-BASED HAIR TEST	<b>Hair Testing</b>	Seven drugs of abuse	<ol style="list-style-type: none"> <li>As early as two days after use.</li> <li>Up to 90 days after use.</li> <li>Depends upon individual metabolism and drug.</li> </ol>	Collection site.	On-site, typically 10 minutes.	As soon as three days	Longest window of detection.